

Product Datasheet

HSP90 alpha (phospho-Thr5/7) antibody (FITC) (orb8337)

Description	HSP90 alpha (phospho-Thr5/7) antibody (FITC)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	FITC
Tested Applications	IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human HSP90 around the phosphorylation site of Thr5 + Thr7
Target	Phospho-HSP90 alpha (Thr5 + Thr7)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P07900
Expiration Date	12 months from date of receipt.