

Product Datasheet

TRPA antibody (orb833790)

Description	TRPA antibody
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Pyrococcus furiosus (strain ATCC 43587 / DSM 3638 / JCM 8422 / Vc1) trpA protein
Target	trpA
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8U094
Expiration Date	12 months from date of receipt.