

Product Datasheet

Trappc2l antibody (orb833956)

Description	Trappc2l antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Mus musculus (Mouse) Trappc2l protein
Target	Trappc2l
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9JME7
Expiration Date	12 months from date of receipt.