

Product Datasheet

UGPE antibody (orb834256)

Description	UGPE antibody
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Escherichia coli (strain UT189 / UPEC) ugpE protein
Target	ugpE
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q1R5H7
Expiration Date	12 months from date of receipt.