

Product Datasheet

UVRC antibody (orb835203)

Description	UVRC antibody
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73) uvrC protein
Target	uvrC
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	B7LPF5
Expiration Date	12 months from date of receipt.