

Product Datasheet

Protein gp55 55 antibody (orb835493)

Description	Protein gp55 55 antibody
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Escherichia phage N15 (Bacteriophage N15) 55 protein
Target	55
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O64363
Expiration Date	12 months from date of receipt.