

Product Datasheet

Y05F antibody (orb835522)

Description	Y05F antibody
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Enterobacteria phage T4 (Bacteriophage T4) y05F protein
Target	y05F
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P39261
Expiration Date	12 months from date of receipt.