

Product Datasheet

Single-stranded DNA-binding protein antibody (orb835526)

Description	Single-stranded DNA-binding protein antibody
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Enterobacteria phage RB8 (Bacteriophage RB8) 32 protein
Target	32
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O21952
Expiration Date	12 months from date of receipt.