

Product Datasheet

Putative DNA transfer protein p33 33 antibody (orb835599)

Description	Putative DNA transfer protein p33 33 antibody
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Acyrthosiphon pisum secondary endosymbiont phage 1 (Bacteriophage APSE-1) 33 protein
Target	33
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9T1R5
Expiration Date	12 months from date of receipt.