

Product Datasheet

YPFJ antibody (orb836874)

Description	YPFJ antibody
Species/Host	Rabbit
Reactivity	E. coli
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC) ypfj protein
Target	ypfj
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P64430
Expiration Date	12 months from date of receipt.