

Product Datasheet

PSAC antibody (orb837136)

Description	PSAC antibody
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Gossypium barbadense (Sea-island cotton) (Egyptian cotton) psaC protein
Target	psaC
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	A0ZZ86
Expiration Date	12 months from date of receipt.