

Product Datasheet

RPL41 antibody (orb837187)

Description	RPL41 antibody
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Gossypium hirsutum (Upland cotton) (Gossypium mexicanum) RPL41 protein
Target	RPL41
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P62122
Expiration Date	12 months from date of receipt.