

Product Datasheet

Gibberellin 20 oxidase 1 antibody (orb839731)

Description	Gibberellin 20 oxidase 1 antibody
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Oryza sativa subsp. japonica (Rice) GA20OX1 protein
Target	GA20OX1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P93771
Expiration Date	12 months from date of receipt.