

Product Datasheet

NSP4 antibody (orb841265)

Description	NSP4 antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Rice gall dwarf virus (RGDV) Non-structural protein 4 protein
Target	Non-structural protein 4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	A4PBP7
Expiration Date	12 months from date of receipt.