

Product Datasheet

NRAMP2 antibody (orb841286)

Description	NRAMP2 antibody
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Oryza sativa subsp. japonica (Rice) NRAMP2 protein
Target	NRAMP2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q10Q65
Expiration Date	12 months from date of receipt.