

Product Datasheet

RPA2A antibody (orb841914)

Description	RPA2A antibody
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Oryza sativa subsp. japonica (Rice) RPA2A protein
Target	RPA2A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6K9U2
Expiration Date	12 months from date of receipt.