

Product Datasheet

NEF antibody (orb842942)

Description	NEF antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human immunodeficiency virus type 1 group M subtype B (isolate HXB2) (HIV-1) nef protein
Target	nef
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P04601
Expiration Date	12 months from date of receipt.