

Product Datasheet

GGAG antibody (orb843260)

Description	GGAG antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Feline sarcoma virus (strain Gardner-Arnstein) (Ga-FeSV) (Gardner-Arnstein feline leukemia oncovirus B) Glyco-Gag protein
Target	Glyco-Gag
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P0DOH1
Expiration Date	12 months from date of receipt.