

Product Datasheet

GAG antibody (orb843271)

Description	GAG antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Feline immunodeficiency virus (isolate TM2) (FIV) gag protein
Target	gag
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P31821
Expiration Date	12 months from date of receipt.