

Product Datasheet

gag-pro-pol antibody (orb843432)

Description	gag-pro-pol antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human T-cell leukemia virus 3 (strain Pyl43) (HTLV-3) gag-pro-pol protein
Target	gag-pro-pol
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q4U0X6
Expiration Date	12 months from date of receipt.