

Product Datasheet

gag-pro-pol antibody (orb843477)

Description	gag-pro-pol antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human T-cell leukemia virus 1 (isolate Caribbea HS-35 subtype A) (HTLV-1) gag-pro-pol protein
Target	gag-pro-pol
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P14078
Expiration Date	12 months from date of receipt.