

Product Datasheet

V-SEA antibody (orb843646)

Description	V-SEA antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Avian erythroblastosis virus (strain S13) V-SEA protein
Target	V-SEA
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P23049
Expiration Date	12 months from date of receipt.