

Product Datasheet

V-SIS antibody (orb843748)

Description	V-SIS antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Woolly monkey sarcoma virus (WMSV) (Smian sarcoma-associated virus) V-SIS protein
Target	V-SIS
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P01128
Expiration Date	12 months from date of receipt.