

Product Datasheet

gag-pro antibody (orb843797)

Description	gag-pro antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Mason-Pfizer monkey virus (MPMV) (Simian Mason-Pfizer virus) gag-pro protein
Target	gag-pro
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P07570
Expiration Date	12 months from date of receipt.