

Product Datasheet

ZC3HAV1 antibody (orb843966)

Description	ZC3HAV1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Homo sapiens (Human) ZC3HAV1 protein
Target	ZC3HAV1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q7Z2W4
Expiration Date	12 months from date of receipt.