

Product Datasheet

Cyp2d4 antibody (orb845325)

Description	Cyp2d4 antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Rattus norvegicus (Rat) Cyp2d4 protein
Target	Cyp2d4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q64680
Expiration Date	12 months from date of receipt.