

Product Datasheet

CZF1 antibody (orb845430)

Description	CZF1 antibody
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Candida albicans (strain SC5314 / ATCC MYA-2876) (Yeast) CZF1 protein
Target	CZF1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q5A0W9
Expiration Date	12 months from date of receipt.