

Product Datasheet

EXO5 antibody (orb845874)

Description	EXO5 antibody
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Saccharomyces cerevisiae (strain Lalvin EC1118 / Prise de mousse) (Baker's yeast) EXO5 protein
Target	EXO5
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	D3UEQ8
Expiration Date	12 months from date of receipt.