

Product Datasheet

Csnk1g3 antibody (orb847449)

Description	Csnk1g3 antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Rattus norvegicus (Rat) Csnk1g3 protein
Target	Csnk1g3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q62763
Expiration Date	12 months from date of receipt.