

Product Datasheet

Map1lc3b antibody (orb847714)

Description	Map1lc3b antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Rattus norvegicus (Rat) Map1lc3b protein
Target	Map1lc3b
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q62625
Expiration Date	12 months from date of receipt.