

Product Datasheet

NUC1 antibody (orb848357)

Description	NUC1 antibody
Species/Host	Rabbit
Reactivity	Yeast
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Saccharomyces cerevisiae (Baker's yeast) NUC1 protein (1-329aa)
Target	NUC1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P08466
Expiration Date	12 months from date of receipt.