

Product Datasheet

rntA antibody (orb849393)

Description	rntA antibody
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Aspergillus oryzae (strain ATCC 42149 / RIB 40) (Yellow koji mold) rntA protein
Target	rntA
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P00651
Expiration Date	12 months from date of receipt.