

Product Datasheet

RVB2 antibody (orb849865)

Description	RVB2 antibody
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Candida albicans (strain SC5314 / ATCC MYA-2876) (Yeast) RVB2 protein
Target	RVB2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q5AGZ9
Expiration Date	12 months from date of receipt.