



Product Datasheet

Histamine 3 Receptor antibody [Out of stock] (orb85050)

Description	Rabbit polyclonal antibody to Histamine 3 Receptor
Species/Host	Rabbit
Reactivity	Equine, Hamster, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Synthetic 16 amino acid peptide from 2nd cytoplasmic domain of human HRH3 / Histamine H3 Receptor. Percent identity with other species by BLAST analysis: Human, Gorilla, Gibbon, Monkey, Mouse, Rat, Hamster, Elephant, Bovine, Horse, Guinea pig (100%); Panda, Bat, Dog (94%); Opossum (81%).
Target	HRH3
Preservatives	PBS, 0.1% Sodium Azide
Concentration	1 mg/ml
Storage	Aliquot and store undiluted at -20°C or below for up to 1 year. Can be stored undiluted at 4°C for up to 1 month. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Entrez	11255
Dilution Range	IHC-P: 7-15 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...



Immunohistochemical staining of human Br...