

Product Datasheet

Histamine H4 Receptor antibody [Out of stock] (orb85089)

Description Rabbit polyclonal antibody to Histamine H4 Receptor

Species/Host Rabbit

Reactivity Human, Monkey

Conjugation Unconjugated

Tested Applications IHC-P

Immunogen Synthetic 18 amino acid peptide from 1st cytoplasmic domain of human HRH4 / Histamine H4 Receptor. Percent identity with other species by BLAST analysis: Human, Monkey, Marmoset (100%); Bovine, Horse, Rabbit (94%); Gibbon, Mouse, Panda, Dog (89%); Rat, Opossum, Platypus (83%).

Target HRH4

Preservatives PBS, 0.1% Sodium Azide

Concentration 1 mg/ml

Storage Aliquot and store undiluted at -20°C or below for up to 1 year. Can be stored undiluted at 4°C for up to 1 month. Avoid freeze-thaw cycles.


Note For research use only

Clonality Polyclonal


Entrez **59340**

Dilution Range IHC-P: 12-24 µg/ml

Expiration Date 12 months from date of receipt.



Immunohistochemical staining of paraffin...



Immunohistochemical staining of human Ov...