

Product Datasheet

Acetyl-Histone H4 (Lys8) antibody (FITC) (orb8510)

Description	Acetyl-Histone H4 (Lys8) antibody (FITC)
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	FITC
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated Synthesised acetylpeptide derived from human Histone H4 around the acetylation site of Lys8
Target	Histone H4 (Acetyl K8)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P62805
Dilution Range	Flow-Cyt=1ug/test,IF=1:100-500,ICC=1:100
Expiration Date	12 months from date of receipt.