


Product Datasheet

Dopamine Receptor D3 antibody [Out of stock] (orb85432)

Description	Rabbit polyclonal antibody to Dopamine Receptor D3
Species/Host	Rabbit
Reactivity	Human, Monkey
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Synthetic 16 amino acid peptide from 3rd cytoplasmic domain of human DRD3 / Dopamine Receptor D3. Percent identity with other species by BLAST analysis: Human, Chimpanzee, Gorilla, Monkey, Marmoset (100%); Gibbon, Mouse, Dog, Bovine, Bat, Hamster, Elephant, Panda, Horse, Rabbit (94%); Rat, Pig, Opossum (88%).
Target	DRD3
Preservatives	PBS, 0.1% Sodium Azide
Concentration	1 mg/ml
Storage	Aliquot and store undiluted at -20°C or below for up to 1 year. Can be stored undiluted at 4°C for up to 1 month. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Entrez	1814
Dilution Range	IHC-P: 20 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...