

Product Datasheet

SPBC13E7.08c antibody (orb856400)

Description	SPBC13E7.08c antibody
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast) SPBC13E7.08c protein
Target	SPBC13E7.08c
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9P6R2
Expiration Date	12 months from date of receipt.