


Product Datasheet

Prostaglandin E receptor 4 antibody [Out of stock] (orb85667)

Description	Rabbit polyclonal antibody to Prostaglandin E receptor 4
Species/Host	Rabbit
Reactivity	Human, Monkey
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Synthetic 17 amino acid peptide from 3rd extracellular domain of human PTGER4 / EP4. Percent identity with other species by BLAST analysis: Human, Chimpanzee, Gorilla, Gibbon, Monkey, Marmoset (100%); Guinea pig (88%); Elephant, Bat, Rabbit, Horse (82%).
Target	PTGER4
Preservatives	PBS, 0.1% Sodium Azide
Concentration	0.9 mg/ml
Storage	Aliquot and store undiluted at -20°C or below for up to 1 year. Can be stored undiluted at 4°C for up to 1 month. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Entrez	5734
Dilution Range	IHC-P: 9-18 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...