

Product Datasheet

TOM5 antibody (orb857655)

Description	TOM5 antibody
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast) tom5 protein
Target	tom5
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	G2TRP7
Expiration Date	12 months from date of receipt.