

Product Datasheet

MAOM antibody (orb858502)

Description	MAOM antibody
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Solanum tuberosum (Potato) NAD-dependent malic enzyme 62 kDa isoform, mitochondrial protein
Target	NAD-dependent malic enzyme 62 kDa isoform, mitochondrial
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P37221
Expiration Date	12 months from date of receipt.