

Product Datasheet

VDAC1 antibody (orb858743)

Description	VDAC1 antibody
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Solanum tuberosum (Potato) Mitochondrial outer membrane protein porin of 34 kDa protein
Target	Mitochondrial outer membrane protein porin of 34 kDa
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P42055
Expiration Date	12 months from date of receipt.