

Product Datasheet

si:dkey-97o5.1 antibody (orb859041)

Description	si:dkey-97o5.1 antibody
Species/Host	Rabbit
Reactivity	Fish
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Danio rerio (Zebrafish) (Brachydanio rerio) si:dkey-97o5.1 protein
Target	si:dkey-97o5.1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q5TYP4
Expiration Date	12 months from date of receipt.