

Product Datasheet

cfap157 antibody (orb859135)

Description	cfap157 antibody
Species/Host	Rabbit
Reactivity	Fish
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Danio rerio (Zebrafish) (Brachydanio rerio) cfap157 protein
Target	cfap157
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	A2BDR7
Expiration Date	12 months from date of receipt.