

Product Datasheet

ctdnep1a antibody (orb859295)

Description	ctdnep1a antibody
Species/Host	Rabbit
Reactivity	Fish
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Danio rerio (Zebrafish) (Brachydanio rerio) ctdnep1a protein
Target	ctdnep1a
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q5U395
Expiration Date	12 months from date of receipt.