

Product Datasheet

eef1akmt1 antibody (orb859508)

Description	Rabbit polyclonal antibody to eef1akmt1. All of our unique antibodies are supplied with the following...
Species/Host	Rabbit
Reactivity	Fish
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Danio rerio (Zebrafish) (Brachydanio rerio) eef1akmt1 protein
Target	eef1akmt1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6NYP8
Expiration Date	12 months from date of receipt.