

Product Datasheet

E4F1 antibody (orb859544)

Description	Rabbit polyclonal antibody to E4F1. All of our unique antibodies are supplied with the following: 200...
Species/Host	Rabbit
Reactivity	Fish
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Danio rerio (Zebrafish) (Brachydanio rerio) e4f1 protein
Target	e4f1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q4V8R6
Expiration Date	12 months from date of receipt.