

Product Datasheet

foxi1e antibody (orb859720)

Description	Rabbit polyclonal antibody to foxi1e. All of our unique antibodies are supplied with the following: 2...
Species/Host	Rabbit
Reactivity	Frog
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Xenopus tropicalis (Western clawed frog) (Silurana tropicalis) foxi1e protein
Target	foxi1e
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6P8A3
Expiration Date	12 months from date of receipt.