

Product Datasheet

gapdh-2 antibody (orb859746)

Description	Rabbit polyclonal antibody to gapdh-2. All of our unique antibodies are supplied with the following: ...
Species/Host	Rabbit
Reactivity	Fish
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Danio rerio (Zebrafish) (Brachydanio rerio) gapdh-2 protein
Target	gapdh-2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q5MJ86
Expiration Date	12 months from date of receipt.