

Product Datasheet

gtpbp10 antibody (orb859830)

Description	Rabbit polyclonal antibody to gtpbp10. All of our unique antibodies are supplied with the following: ...
Species/Host	Rabbit
Reactivity	Fish
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Danio rerio (Zebrafish) (Brachydanio rerio) gtpbp10 protein
Target	gtpbp10
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6DHF7
Expiration Date	12 months from date of receipt.